



TO PAN


Dental Panoramic Phantom

A test object that allows the complex focal trough of the panoramic x-ray unit to be quantified, as well as the capacity for spatial resolution testing using line-pair test patterns (available separately).

This allows full assessment of parameters affecting image quality.
Can be used in film and digital systems.

TO PAN Enables:

- Quantification of the focal trough and image layer
- Accurate determination of the bite offset
- Quantification of the symmetry of the scan
- Measurement of the incident beams angles
- Simple visualisation of faults
- Spatial resolution testing



For more information visit us at www.leadstestobjects.com



TO PAN

Dental Panoramic Phantom



fig. 1 Radiograph of TO PAN

