US-10 US-10A US-10B US-10C US-10D


Female Pelvic Ultrasound Phantom

Practice both transvaginal and transabdominal ultrasound procedures in one

Newly added Follicular phase unit support thorough training in Obs/gyn ultrasound


Product Supervision

Charlotte Henningsen, MS, RT(R), RDMS, RVT, FSDMS, FAIUM
Chair & Professor - Sonography Department
Adventist University of Health Sciences





transvaginal



transabdominal

FEATURES

1 | Compatible with actual ultrasound transducer
2 | Six types of interchangeable phantom inserts with different pathologies
3 | Realistic view of female external genital
4 | Excellent ultrasound image quality

SKILLS


| Handling and manipulation of transvaginal and transabdominal transducers
| Interpretation of sonogram
| Visualization and localization of anatomies and pathologies

ANATOMY

Common for all units
| uterus
| urinary bladder
| pelvic bone
| ovaries
| follicles
| vagina
| rectum





Pregnancy units
| gestation sac
| yolk sac
| embryo
| placenta
| follicles
| umbilical cord (10-week)

longitudinal view of the pregnancy (10-week) unit on transvaginal scan



For set variations, please refer the website:
https://www.kyotokagaku.com/en/products_introduction/us-10/

KEY FEATURES

Six variations of training units	Anatomy, findings and pathology
Normal pregnancy (7-week) unit	
Normal pregnancy (10-week) unit	
Ectopic pregnancy unit •ectopic pregnancy •bleeding at Douglas cavity	
Pathological unit •endometrial cancer •dermoid cyst •uterine fibroid •bleeding at Douglas cavity	
IUD unit	
NEW Follicular phase unit	

SPECIFICATIONS

Torso size: W34 x D33 x H24 cm
W13.5 x D13 x H9.5 in.

weight: 6.9kg / 15.2 lbs

MATERIALS

Soft resin, Polyurethane elastomer
Latex free

* Ultrasound scanner is not included.

DESCRIPTIONS

REPLACEMENT PARTS

41909-030	ultrasound pathology unit
41909-040	ultrasound ectopic pregnancy unit
41909-050	ultrasound first trimester pregnancy unit (7W)
41909-060	skin for US-10 / 10A / 10B
41909-070	ultrasound first trimester pregnancy unit (10W)
41909-080	IUD unit
41909-090	follicular phase unit
41909-300	lower torso with accessories